

GUARNIFLON S.p.A. - PATI Division

Product code : OTQHST9900 - 200 mic

Direction : TD

Preload : 0,05 MPa


Rate of grip separation : 500 mm/min

Initial grip separation : 50,00 mm

Risultati del test:

Sample no.	Thickness μm	S_0 mm^2	E_t 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
1	196	1,96	884	21,6	56,6	370
2	196	1,96	907	21,7	62,4	400
3	195	1,95	905	22,6	55,4	360
4	200	2,00	937	22,2	52,4	350

Grafico della serie:



Statistiche:

Serie	Thickness μm	S_0 mm^2	E_t 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
n = 4						
\bar{x}	196,8	1,97	908	22,0	56,7	370