

GUARNIFLON S.p.A. - PATI Division

Product code : OTQHST9900 - 250 mic

Direction : TD

Preload : 0,05 MPa


Rate of grip separation : 500 mm/min

Initial grip separation : 50,00 mm

Risultati del test:

Sample no.	Thickness μm	S_0 mm^2	E_t 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
1	257	2,57	719	23,8	50,2	380
2	255	2,55	845	23,3	51,2	380
3	255	2,55	849	23,2	49,8	380
4	261	2,61	894	23,9	51,5	380
5	259	2,59	865	23,4	49,9	380

Grafico della serie:



Statistiche:

Serie	Thickness μm	S_0 mm^2	E_t 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
n = 5						
\bar{x}	257,4	2,57	834	23,5	50,5	380